CYTOPLAST®

RTM Collagen

INSTRUCTIONS FOR USE

EN

English

DE

Deutsch

FR

Français

ES

Español

IT

Italiano

PT

Português

SV

Svenska

DA

Dansk

NL

Nederlands

CYTOPLAST[®]

RTM Collagen - Resorbable Collagen Membrane

INSTRUCTIONS FOR USE

Introduction: The user of Cytoplast[®] products has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Osteogenics Biomedical disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use or installation of Cytoplast[®] products.

Intended Use: Cytoplast[®] RTM Collagen is a bioresorbable, implantable collagen material that is intended for use in dental surgery procedures as a material for placement in the area of dental implant, bone defect or ridge reconstruction to aid in wound healing post dental surgery.

Description: Cytoplast[®] RTM Collagen is a white, nonfriable membrane matrix engineered from highly purified type I collagen fibers derived from bovine Achilles tendon. Cytoplast[®] RTM Collagen is resorbable which eliminates the need for a second surgical procedure that is normally required to remove a non-resorbable membrane. Cytoplast[®] RTM Collagen is supplied sterile, non-pyrogenic and for single use only.

Administration: Cytoplast[®] RTM Collagen is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments.

The periodontal defect or bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). The clinician should perform thorough debridement of the defect. Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect. As much tissue as possible should be preserved to allow for primary closure of the wound and correct positioning of the flaps without tension.

Cytoplast[®] RTM Collagen can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to the final placement. Then, the device may be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors. The membrane should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. Primary closure must be achieved and maintained to avoid accelerated resorption due to membrane exposure.

Post-operative Procedures: Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent. The patient should refrain from brushing the treated area for two days following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated. The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

Cytoplast[®] RTM Collagen should be completely resorbed 26 to 38 weeks following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

Contraindications: Cytoplast[®] RTM Collagen is contraindicated in patients who have: acute infections or contaminated wound in the oral cavity, known allergy to collagen of animal origin or other bovine-derived products, clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment, will prevent safe implantation or likely healing.

Warning: Clinicians should use care in screening their patients for any known allergies to collagen or bovine-derived products. Hypersensitivity reactions have been noted with the use of other products containing bovine collagen; therefore, the possibility exists of developing a local sensitivity response to Cytoplast® RTM Collagen.

Precautions: As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, indicating a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. Nor has it been evaluated in pregnant women, children and/or in patients with conditions involving extremely severe defects with little peridontium or bone. Cytoplast® RTM Collagen cannot be re-sterilized. Open, unused membranes must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used. Do not use if the product sterilization barrier or its packaging is compromised.

Adverse Reactions: Possible complications that can occur with any dental surgery include infection, swelling of the intraoral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing or dehiscence, root resorption or ankylosis, loss of crestal bone height, pain, or complications associated with the use anesthesia. Minor discomfort may occur for a few days.

Safety: The product is manufactured from bovine Achilles tendon, which is classified as tissues with no detected infectivity for Bovine Spongiform Encephalopathy, BSE (World Health Organization Guidelines). The bovine tendon is known to be one of the richest sources of type I collagen that is commercially available.

The manufacturing process for the product meets European and International Standards for animal tissue sourcing, handling and inactivation of Spongiform Encephalopathy (SE) pathogens. This process involves a treatment with sodium hydroxide that is a recognized method of inactivation of SE pathogens.

A viral inactivation study for the product's manufacturing process was conducted by an independent laboratory. In this study key manufacturing steps were evaluated for their ability to inactivate the following viral strains: Bovine Viral Diarrhea (enveloped virus) and Porcine Parvoviridae (non-enveloped virus). The study results showed that each of the manufacturing steps evaluated, including the sodium hydroxide treatment, is effective in inactivating these viruses.

Storage: The product should be stored at room temperature. Avoid excessive heat and humidity.

How Supplied: One (1) membrane per package; Two (2) packages per box.

Available Sizes:	RTM1520	15 mm x 20 mm
	RTM2030	20 mm x 30 mm
	RTM3040	30 mm x 40 mm

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or dentist.

Labeling Symbols Symbols may be used on package labeling for easy identification.

Manufacturer

Catalog Number

Use by

Do not use if the product sterilization barrier or its packaging is compromised.

Do Not Reuse

Temperature limitations – store at room temp. 15-30° C (59-86°F)

Attention, see instructions for use

Method of Sterilization Using Irradiation

Lot Number

Authorized Representative in the European Community

Manufacturer: Collagen Matrix, Inc.
15 Thornton Road
Oakland, NJ 07436 USA

Distributed by: Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424 USA
www.osteogenics.com

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Fabrikant: Collagen Matrix, Inc.
Gedistribueerd door: Osteogenics Biomedical, Inc.
MDSS GmbH