ZCORE™

porcine xenograft particulate

INSTRUCTIONS FOR USE

EN

English

ES

Español

SV

Svenska

DE

Deutsch

IT

Italiano

DA

Dansk

FR


Français

PT

Português

NL

Nederlands

EN

ZCORE™ Porcine Xenograft Particulate

A sterile, biocompatible anorganic porous bone mineral for use in periodontal, oral and maxillofacial surgery.

DESCRIPTION

ZCORE™ is a porous bone mineral matrix consisting predominantly of calcium phosphate. It is produced by removal of the organic components from porcine cancellous bone. ZCORE™ is sterilized by gamma irradiation. The product is available in granular form, nonpyrogenic and for single use only.

PROPERTIES/ACTIONS

The anorganic bone matrix of ZCORE™ has interconnecting macro and microscopic porous structure that supports the formation and ingrowth of new bone at the implantation site. The use of ZCORE™ may be considered when autogenous bone is not indicated, or insufficient in quantity to fulfill the needs of the proposed surgical procedure.

INDICATIONS AND USAGE

ZCORE™ is indicated for:

- Augmentation or reconstructive treatment of the alveolar ridge.
- Filling of infrabony periodontal defects.
- Filling of defects after root resection, apicoectomy, and cystectomy.
- Filling of extraction sockets to enhance preservation of the alveolar ridge.
- Elevation of the maxillary sinus floor.
- Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR).
- Filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR).

INSTRUCTIONS FOR USE

- After exposure of the bony defect with a mucoperiosteal flap, all granulation tissue must be carefully removed.
- Mix ZCORE™ with autogenous bone, osseous coagulum, patient's blood or sterile normal saline. If large maxillofacial defects are present, ZCORE™ should be mixed with autogenous bone in a ratio of approximately 1:1.
- In order to assure the formation of new bone, ZCORE™ should only be placed in direct contact with well vascularized bone. Cortical bone should be mechanically perforated.
- Loosely pack ZCORE™ into osseous defect using a sterile instrument. Use of excessive force

- will result in compression of the particles and loss of trabecular architecture.
- Overfilling of the defects should be avoided.
- The mucoperiosteal flaps should be sutured to achieve primary closure, if possible. A surgical dressing may be placed over the wound for one to two weeks.
- If primary wound closure cannot be achieved completely, further immobilization of the flap (e.g., by incision through the periosteum) should be performed and/or a bioabsorbable membrane (e.g. Resorbable Collagen Membrane) should be placed over the bone graft site.
- Post-operative procedures: Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for 2 weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent. Consult the manufacturer's Instructions for Use for side effects and other known risks associated with the use of chlorhexidine gluconate oral rinses.

CONTRAINDICATIONS

Contraindications customary to the use of bone grafts should be observed. *ZCORE™* should not be used in patients with:

- Acute or chronic infection (osteomyelitis) at the surgical site
- Metabolic diseases (diabetes, hyperparathyroidism, osteomalacia)
- Severe renal dysfunction, severe liver disease
- High dose corticosteroid therapy
- Vascular impairment at the implant site
- Osteoporosis
- Known allergy or hypersensitivity to porcine-derived implantable materials

WARNINGS

The device should be secured to prevent motion and migration, use in areas where the graft can be adequately contained.

Do not use if package is opened or damaged or if expiration date has been exceeded.

ZCORE™ cannot be re-sterilized or re-used. Open, unused product must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

PRECAUTIONS

In order to facilitate the formation of new bone *ZCORE™* should only be implanted in direct contact with a well vascularized bony tissue. Drilling may be recommended to facilitate bleeding from cortical bone.

In larger defects a mixture of autogenous bone or bone marrow may improve the formation of new bone.

Implantology

Generally, in augmented areas, the placement of titanium fixtures should take place once the bone has sufficient strength and integrity for dental implant placement, which is typically greater than 6 months after implantation of a bone graft material. For sinus floor elevation, typically 9-12 months should be allowed after implantation of bone graft material before placement of the titanium fixtures. Xrays should be taken to confirm the bone integrity prior to dental implant placement.

Periodontology

The filling of periodontal defects with ZCORE™ requires (along with plaque control) the successful local treatment of the periodontal lesion (e.g. root planing, debridement of granular tissue) prior to implantation.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist.

ADVERSE REACTIONS

Possible complications that may occur with any dental surgery include swelling at the surgical site, flap sloughing, bleeding, local inflammation, bone loss, infection or pain.

HOW SUPPLIED

ZCORE™ is supplied sterile, non-pyrogenic, and for single use only.

Catalogue Number	Particle Size	Volume
ZS050	0.25 – 1.0 mm	0.5 cc
ZS100		1.0 cc
ZS200		2.0 cc
ZS400		4.0 cc
ZL100	1.0 – 2.0 mm	1.0 cc
ZL200		2.0 cc

STORAGE

Store at room temperature (15°C/59°F - 30°C/86°F).

LABELING SYMBOLS

Symbols may be used on some international package labeling for easy identification.

Caution

Use By

Do Not Reuse

LOT

Lot Number

STERILE R

Sterilized Using Irradiation

REF

Catalog Number

R_x only

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist

Manufacturer

Do not use if the product sterilization barrier or its packaging is compromised.

Temperature limitation

EC REP

Authorized Representative in the European Community

Manufacturer:
Collagen Matrix, Inc.
15 Thornton Road
Oakland, New Jersey 07436 USA

Distributed by:
Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424 USA
www.osteogenics.com
+1-806-796-1923

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MDSS GmbH