

# Master-Pin-Control

*revolutionary hybrid pin system*



*The Master-Pin-Control Bone Management® system is used for the fixation of membranes (absorbable and non-absorbable) in order to avoid micro-mobility of the graft. The pins have an extremely sharp tip that allows precise placement into cortical bone. Mini-threads on the pins make them a hybrid of a screw and pin. The threads on the pins increase the surface area of the shaft, resulting in pin stability, while also making removal of the pins possible with the included screwdriver.*

*Meisinger* since 1888



### Master-Pin-Control

| BMP00

**\$1,495**

(34) Pins

Master-Pin-Tray

Screw Driver For Pin Removal

Fixation Holder

Initial Bur

Twist Drills: (2) 0.6 mm twist drills, (2) 0.8 mm twist drills



### Master-Pin-Control Basic

| BMPBA

**\$1,295**

(10) Pins

Master-Pin-Tray

Screw Driver For Pin Removal

Fixation Holder

Initial Bur

Twist Drills: (2) 0.6 mm twist drills, (2) 0.8 mm twist drills



### Master-Pin Longer Screw Driver

*Designed to make pin removal easier in hard to reach areas*

| MP15 (1 per box)

**\$157**



### Replacement Pins

| MP10 (10 per box)

**\$180**

| MP50 (50 per box)

**\$734**



### Decortication Bur

*1.2 mm diameter x 4.0 mm long decortication bur with drill stop*

| 203S-012-RA (2 per box)

**\$90**