

Safescraper® Twist - Curve Version

Versatile cortical bone collector

Holds up to 2.5 cc at a time

3987 (3 sterile scrapers per package)



Applications

- Extraction defects
- Periodontal defects
- Sinus lift procedures
- Ridge augmentation

Harvesting Sites

- Oblique external line w/ tunnel
- Ramus
- Mandibular symphysis
- Sinus window
- Lingual bone
- Zygomatic area
- Nasal spine
- Palate
- Small areas near the defect

Features & Benefits

Ergonomic design

Cortical bone harvesting is easily achieved from intraoral sites with a minimally invasive approach

2.5 cc collection chamber

Large amounts of bone may be collected at once

Bone is collected with coagulated blood

Graft has high biological plasticity, making it easy to handle and mold


Superior harvesting method

The manual harvesting technique allows graft to retain cell viability that can be compromised with other harvesting techniques that mill, grind, or potentially overheat bone

Safe

The disposable scraper is sterile and allows clinicians to harvest autogenous bone, which eliminates any chance of disease transmission

A 160° blade allows clinicians to collect bone from any bony surface.



The Safescraper® Twist's transparent chamber holds up to 2.5 cc of bone, which can be used alone or mixed in combination with other graft materials.



"This unit works really well and has nice contours to use in difficult harvesting sites."

Tom Faerber, DMD; Oral and Maxillofacial Surgeon