

Smartscraper

Cortical bone collector and syringe in one



Holds up to 0.3 cc at a time


4890 (3 sterile scrapers per package)

Applications

- Extraction defects
- Periodontal defects
- Sinus lift procedures
- Ridge augmentation

Harvesting Sites

- Oblique external line w/ tunnel
- Ramus
- Mandibular symphysis
- Sinus window
- Lingual bone
- Zygomatic area
- Nasal spine
- Palate
- Small areas near the defect



The Smartscraper is opened with a simple movement. The syringe, in which the bone particulate has been collected, can then be used to place graft directly into areas with limited access.

